## Binokularmikroskope

MBL2000-Serie

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## Binocular Microscope

MBL2000 Series

Operating Manual

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## Index

### ENGLISH

|                                         |    |
|-----------------------------------------|----|
| Introduction                            | 31 |
| Models Overview                         | 32 |
| 1. Description MBL2000                  | 34 |
| 2. Assembly of the microscope           | 36 |
| 3. Operating the microscope             | 37 |
| 3.1 Eye distance and dioptre adjustment | 38 |
| 3.2 Coarse adjustment                   | 38 |
| 3.3 Aperture diaphragm                  | 39 |
| 3.4 Oil immersion lenses                | 39 |
| 3.5 Cover glass and object holder       | 40 |
| 4. Technical Data                       | 41 |
| 5. Definitions                          | 42 |
| 6. Accessories                          | 43 |
| 6.1 Object lenses                       | 43 |
| 6.2 Eyepieces                           | 43 |
| 7. Phase contrast device                | 44 |
| 8. Dark field condenser                 | 50 |
| 9. Photo- and video mounting            | 51 |
| 10. Maintenance                         | 52 |
| Recovery and recycling                  | 53 |
| Warranty conditions                     | 54 |

## Introduction

Congratulations on your purchase of a new microscope!

The microscope MBL2000 series has been developed for universal applications in research, medicine and industry. The high manufacturing quality as can be seen, for example, with the metal casing, as well as the high-quality optical components ensure a long product life. Thanks to the wide range of accessories available, these microscopes can be used for many applications.

All microscopes of the MBL2000 series consist of a stable metal stand and have been equipped with a 45° inclined optical head that can be rotated by 360°. The optical system consists of a 10x flat-field eyepiece and interchangeable achromatic object lenses. 20x eyepiece or flat-field achromatic object lenses for especially demanding applications are also available. The illumination consists of a 20 or 30 Watt, 6V halogen bulb, the bright field Abbe condenser and the iris diaphragm. A phase contrast device, a dark field feature or an optical head with phototube is also available.

A. KRÜSS would like to point out that this instruction manual contains important safety and maintenance information and should therefore be made available to all users. A.KRÜSS does not accept any liability for the improper use of the microscope.

### Important handling instructions

- The microscope is a precision instrument. Handle it with the appropriate care.
- Protect you microscope from dust, moisture, vibration and especially from direct exposure to sunlight.
- Never turn the opposing fine and coarse adjustment knobs both at the same time.
- Remove the power plug from the electrical socket before replacing the bulb or fuse.
- Please always take the advices in Chapter 10 | Maintenance into consideration.

## Models Overview


**T** Trinocular microscope / photo tube  
**PL** Flat-field achromatic object lenses  
**PH** Phase contrast device  
**B** Blood test device  
**30W** 30Watt illumination

|                              | Optical features                                                                                                                                 | Equipment                                                                                              | Illumination                                                                | Characteristics                                     | Application                                                |
|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|-----------------------------------------------------|------------------------------------------------------------|
| <b>MBL2000 (basic model)</b> | Binocular-Microscope, 10x flat-field eyepiece, object lenses: 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil                             | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 20 W adjustable, bright field ABBE condenser                            |                                                     | Lab, quality control, universities                         |
| <b>MBL2000-30W</b>           | Binocular-Microscope, 10x flat-field eyepiece, object lenses: 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil                             | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 30 W adjustable, bright field ABBE condenser                            | 30 W illumination                                   | Lab, quality control, universities                         |
| <b>MBL2000-T</b>             | Trinocular-Microscope, 10x flat-field eyepiece, object lenses: 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil                            | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 20 W adjustable, bright field ABBE condenser                            | Photo tube                                          | Lab, quality control, universities                         |
| <b>MBL2000-T-30W</b>         | Trinocular-Microscope, 10x flat-field eyepiece, object lenses: 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 60x/NA 0.85, 100x/NA 1.25 oil               | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 30 W adjustable, bright field ABBE condenser                            | Photo tube<br>30 W illumination                     | Lab, quality control, universities                         |
| <b>MBL2000-PL-PH</b>         | Binocular-Microscope, 10x flat-field eyepiece, object lenses (flat field achrom.): 4x/NA 0.10, PH10x/NA 0.25, PH40x/NA 0.65, PH100x/NA 1.25 oil  | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 20 W adjustable, bright field ABBE condenser phase contrast, dark field | Phase contrast device with dark field               | Sewage treatment plant, lab, quality control, universities |
| <b>MBL2000-T-PL-PH</b>       | Trinocular-Microscope, 10x flat-field eyepiece, object lenses (flat field achrom.): 4x/NA 0.10, PH10x/NA 0.25, PH40x/NA 0.65, PH100x/NA 1.25 oil | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 20 W adjustable, bright field ABBE condenser phase contrast, dark field | Photo tube<br>Phase contrast device with dark field | Sewage treatment plant, lab, quality control, universities |
| <b>MBL2000-PL</b>            | Binocular-Microscope, 10x flat-field eyepiece, object lenses (flat field achrom.): 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil        | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 20 W adjustable, bright field ABBE condenser                            | Flat field achromatic object lenses                 | Lab, quality control, universities                         |

|                         |                                                                                                                                      |                                                                                                        |                                                                                          |                                                                                                                    |                                                                                                 |
|-------------------------|--------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| <b>MBL2000-PL-30W</b>   | Binocular-Microscope, 10x flat field, object lenses (flat field achromatic): 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil  | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 30 W adjustable, bright field ABBE condenser                                         | Flat field achromatic object lenses, 30 W illumination                                                             | Lab, quality control, universities                                                              |
| <b>MBL2000-T-PL</b>     | Trinocular-Microscope, 10x flat field, object lenses (flat field achromatic): 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 20 W adjustable, bright field ABBE condenser                                         | Photo tube<br>Flat field achromatic object lenses                                                                  | Lab, quality control, universities                                                              |
| <b>MBL2000-T-PL-30W</b> | Trinocular-Microscope, 10x flat field, object lenses (flat field achromatic): 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 6 V 30 W adjustable, bright field ABBE condenser                                         | Photo tube,<br>Flat field achromatic object lenses, 30 W illumination                                              | Lab, quality control, universities                                                              |
| <b>MBL2000-B</b>        | Binocular-Microscope, 10x flat field, object lenses: 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil                          | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 15 V 150 W adjustable, cold light source adjustable, dark field condenser for blood test | 150 W cold light source, dark field for blood                                                                      | Blood tests according to Enderlein, alternative practitioner lab, quality control, universities |
| <b>MBL2000-T-B</b>      | Trinocular-Microscope, 10x flat field, object lenses: 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil                         | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 15 V 150 W adjustable, cold light source adjustable, dark field condenser for blood test | Photo tube,<br>150 W Cold light source, dark field for blood, third tube for connecting to photo and video cameras | Blood tests according to Enderlein, alternative practitioner lab, quality control, universities |
| <b>MBL2000-PL-B</b>     | Binocular-Microscope, 10x flat field, object lenses (flat field achromatic): 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil  | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 15 V 150 W adjustable, cold light source adjustable, dark field condenser for blood test | 150 W cold light source, dark field for blood, flat field achromatic object lenses                                 | Blood tests according to Enderlein, alternative practitioner lab, quality control, universities |
| <b>MBL2000-T-PL-B</b>   | Trinocular-Microscope, 10x flat field, object lenses (flat field achromatic): 4x/NA 0.10, 10x/NA 0.25, 40x/NA 0.65, 100x/NA 1.25 oil | XY-stage, coaxial fine and coarse adjustment, iris diaphragm, filter holder, blue filter, green filter | 15 V 150 W adjustable, cold light source adjustable, dark field condenser for blood test | Photo tube,<br>150 W cold light source, dark field for blood, flat field achromatic object lenses                  | Blood tests according to Enderlein, alternative practitioner lab, quality control, universities |

1. Description

1. Description MBL2000

1. Description

## 2. Assembly of the microscope

For packaging purposes, the microscope is delivered in several separate parts. The assembly of the microscope is easy.

1. Remove the stand from the packaging, attach the ocular head and fasten it by using the screw located at the side.

2. Remove the covers and insert the eyepieces into the ocular head.

3. Unpack the four lenses from the plastic boxes and screw them into the lens turret.

The microscope is now ready for use. In the appendix you will find an explanation on how to attach the accessories.

## 3. Operating the microscope

1. Turn on the main switch.
2. Place the sample on the microscope stage.
3. Release the coarse adjustment pre-selection stop (see 3.2).
4. Focus with the 4x or 10x object lens on the sample.
5. Adjust the eye distance and dioptre compensation (see 3.1).
6. Swivel in the desired object lens.
7. Lock the adjusting level for the coarse adjustment pre-selection stop into position (see 3.2).
8. Adjust the light intensity using the brightness control.
9. Fine adjust the focus using the fine adjustment knob.
10. Adjust the aperture diaphragm (see 3.3)

### Note

In case there is no power available, the mirror (available as an accessory) can be used to illuminate the sample. For this purpose, remove the bulb condenser and replace it with the mirror. The base of the mirror is screwed into the socket of the bulb condenser.

### 3.1 Eye distance and dioptre adjustment

1. Take the eyepiece carriage of the right and left eyepiece support with both hands and adjust the proper eye distance
2. The eye distance is correctly adjusted when the field of vision is in focus with the eyes relaxed.
3. Take note of your eye distance shown on the scale.
4. Now turn the longitudinal tube equalizer ring of the right eyepiece support to the value indicated on the eye distance scale.
5. Look through the right eyepiece and focus the object with the fine adjustment knob.
6. Now look through the left eyepiece and adjust the tube equalizer ring without touching the coarse or fine adjustment knob until the sample is in focus.

All transmitted light microscopes have a mechanical tube length of 160 mm.

### 3.2 Coarse adjustment

The coarse adjustment can be adjusted to run smooth or rough by turning the knurled set ring between the microscope stand and the right-hand coarse adjustment knob. This prevents the sliding of the microscope stage and any associated focusing errors.

#### Important note

Never counter-rotate the right and left adjustment knobs at the same time.

A possible contact between the sample and the object lens can be prevented by the pre-selection stop which also makes the coarse focusing easier. After focusing with the coarse adjustment knob, the lever is turned in a clockwise direction up to the stop position. The microscope stage can then no longer be moved beyond this position with the coarse adjustment knob.


Even after the stage is moved to change the object lens, it can only be moved to the pre-selected position with the coarse adjustment knob provided that the fine adjustment knob has not been re-adjusted. The sample is therefore back in focus. The adjustment knob is not affected by the pre-selection stop.

#### Please note!

When adjusting with the coarse adjustment knob, you have to ensure that the pre-selection stop is adjusted to that particular object lens. Otherwise, it is possible that during the simultaneous adjustment and viewing through the eyepiece, the object lens is moved into the sample damaging the object lens as well as the sample.

### 3.3 Aperture diaphragm

The aperture diaphragm integrated in the condenser is used to adjust the depth of sharpness of the image contrast and the resolution power of the microscope lenses. In order to achieve a perfect image quality, the diaphragm opening should generally be open 70-80% of the object lens aperture. You can correctly adjust the aperture diaphragm by removing the eyepiece. You can control the image of the aperture diaphragm by the light path (you might have to dim the light intensity). Theoretically, the aperture diaphragm must be readjusted after each change of the object lens.

### 3.4 Oil immersion lenses, marked with "HI" – homogeneous immersion

To make full use of the numeric aperture of an immersion lens, the object lens and sample – in special cases also the condenser – must be submerged in immersion oil.

### 3. Operating the microscope / 3.5 Cover glass and object holder

Please proceed as follows:

1. First focus on the sample to be examined with a low magnifying lens (e.g. 10x).
2. Use the coarse adjustment pre-selection stop as described in 3.2.
3. Lower the microscope stage by using the coarse adjustment knob.
4. Apply a drop of immersion oil on the spot where the sample is to be examined and which is illuminated by the condenser.
5. Swivel the immersion lens into position. Carefully raise the microscope stage with the coarse adjustment knob as far as it will go (to the pre-selection stop).
6. If needed, refocus with the fine adjustment knob. Make sure that there are no air bubbles developing in the immersion oil.

After use, clean the immersion lens with lens cleaning paper or a suitable cloth using ether or alcohol).

**Important note!**

Make sure that during the swivelling of another object lens, it does not come into contact with the oil, thus becoming contaminated!

### 3.5 Cover glass and object holder

The cover glass is placed over the sample. The lenses provided with an engraving of 0.17 have been corrected for a cover glass density of 0.17 mm, i.e. in case of dry lenses with a NU (numeric aperture) of 0.7 and higher, only a cover glass density of 0.17 mm is to be used as the object lens comes with a cover glass correction.

### 4. Technical data

A special cover glass for blood counting chambers of 0.4 mm thickness may be used for the achromatic 40x NA 0.65 object lens or for object lenses with lower magnifications.

Object holder:

We recommend you use object holders with a thickness between 0.8 and 1.5 mm.

## 4. Technical data

### Optical head

Inclined optical head, symmetrical eye distance adjustment (55 – 75 mm), dioptrre compensation with scale.

### Lens turret

Quadruple

### Achromatic object lenses

4x 0.10 NA 17.04 mm working distance  
10x 0.25 NA 8.05 mm working distance  
40x 0.65 NA 0.32 mm working distance  
100x 1.25 NA 0.13 mm working distance

### Eyepieces

10x, focal length 25 mm

### Stand

Made of metal with coaxial coarse/fine knob with a range of 30 mm each.  
The right coarse adjustment knob comes with a mobility adjustment, the left coarse adjustment knob comes with a quick focusing adjustment.

### Microscope stage

With 0.1 mm increment memory scale, left-right range of movement 74 mm, forwards - backwards 30 mm.

### Low voltage illumination

Integrated, 6 V 20 W (or 30 W) halogen bulb with brightness control.

### Condenser

Double lens ABBE Condenser, NA 1.25, with centring, height adjustment and filter holder.

### CE conformity

The device complies with  
EN 50081-1, IEC 1000-4-2,  
IEC 1000-4-4 und IEC 1000-4-5.

## 5. Definitions

### Working distance

Distance between sample and front lens mount

### Numeric aperture

A mathematical relationship, directly related to the (NA) of the aperture ratio and the resolution power of the object lens.

The NA is the product of the sine of half the aperture angle of an object lens and of the refractive index of the medium through which light enters the object lens.

This is an especially important constant for high-performance object lenses.

### Resolution power

The capacity of the object lens to separate the fine lens structures of the sample. The resolution power of a microscope lens is measured by how two adjacent dots can be separated. (line structures in the object can be seen as a set of points). The resolution power of microscope lens is defined as follows:

$$A = 0.61 \frac{\text{wavelength}}{\text{NA}}$$

The wavelength range of visible light is between 400  $\mu\text{m}$  and 700  $\mu\text{m}$ .

The resolution power increases when the light wavelength decreases. The higher the resolution power of the object lens, the finer the structure of the object can be resolved.

### Depth of sharpness

Distance between the upper and lower limit of the image sharpness which is determined by the optical system.

Structures lying outside this limit appear out of focus and can only be sufficiently seen when using low magnifying lenses. A low depth of sharpness is a special disadvantage in microphotography.

### Field of view number

The number that specifies the diameter of the field in the eyepiece (diaphragm) in mm.

### Diameter of field of view

The actual sample field in mm. The diameter of the field of view of a microscope eyepiece in mm is calculated as follows:

$$\frac{\text{Field of view number of the eyepiece}}{\text{Field of view on the object side in mm}} = \frac{\text{Magnification}}{\text{of the lens}}$$

## 6. Accessories

### 6.1 Object lenses

As a standard, the microscope is equipped with achromatic lenses 4x, 10x, 40x and 100x (oil).

A 20x achromatic object lens is also available for a 200x total magnification, a 60x achromatic lens for a 600x magnification. In order to achieve improved magnification of the perimeters of the field of view, the achromatic lenses can be replaced by flat field achromatic lenses.

The flat field achromatic lenses are also available in 4x, 10x, 20x, 40x and 100x (oil).

Object lens data:

| Vergrößerung | NA   | Arbeitsabstand       |
|--------------|------|----------------------|
| 4x           | 0.10 | 17.04 mm             |
| 10x          | 0.25 | 8.05 mm              |
| 20x          | 0.40 | 2.10 mm              |
| 40x          | 0.65 | 0.32 mm teleskopisch |
| 100x         | 1.25 | 0.13 mm Öl-Immersion |

### 6.2 Eyepieces

As a standard, the microscope is equipped with one or two 10x wide-field flat field eyepieces. This results in a total magnification of 40x, 100x, 400x and 1000x.

When using of one or two 16x eyepieces you can achieve total magnifications of 64x, 160x, 640x and 1600x.

A 10x pointer eyepiece accessory is available for measuring and pointing to samples.

A pointer is directed to the centre is in the image (see fig. 1). By using the scales on the XY stage, it is possible to measure the distances.


Abb. 1

Abb. 2

Another accessory for measuring samples is the micrometer eyepiece (see fig. 2)

The micrometer eyepiece has a 15x magnification which achieves a total magnification of 60x, 150x, 600x and 1500x.

Depending on the object lens in use, the scale length of 10 units corresponds to:

- Object lens 4x = 2500 µm
- Object lens 10x = 1000 µm
- Object lens 40x = 250 µm
- Object lens 100x = 100 µm

Only one pointer or micrometer eyepiece should be used for the binocular microscope as otherwise there maybe a double image of the scale.

## 7. Phase contrast device

The phase contrast device can be adapted to the following microscopes:

MML1500, MBL2000 and MBL2000-T

### MML 2032

The phase contrast device MML 2032 consists of a positive phase contrast lens with a phase lag ring (20x), a centring ring diaphragm and a centring telescope, dark field and bright field.

### MML 2031

The phase contrast device MML 2031 consists of a positive phase contrast lens with a phase lag ring (40x), a centring ring diaphragm and a centring telescope, dark field and bright field.

### MML 2030

The phase contrast device MML2030 consists of a positive phase contrast lens and a centring diaphragm for 10x, 40x and 100x magnifications and a centring telescope.


## Technical data and specifications

The system works according to the positive phase contrast method. The individual phase contrast elements have the following data:

- PH-1** Achromatic phase contrast lens  
10x NA 0.25 with ring diaphragm in a centring mount.
- PH-2** Achromatic phase contrast lens  
40x NA 0.65 with ring diaphragm in a centring mount.
- PH-3** Achromatic phase contrast lens 100x NA 1.25 for oil immersion with diaphragm ring in a centring mount.  
Centring telescope with an adjustable eyepiece.

## Functional principle

When examining colourless, transparent samples in the bright field, the contrast is often too low. In order to avoid colouring or to examining living samples, the phase contrast method can be used. Coloured samples can be examined because the different components have different densities. In case of transparent samples, details in the bright field are not easy to examine as all components have the same density.

Light path when using the phase contrast method

The dark field can be used to obtain sharp images due to border effects such as light scattering and diffraction. If the sample has a crystalline structure or directional characteristics, polarized light can also be used. The phase contrast method is a form of illumination in which one part of the light beam is superimposed on the rest of the light beam to produce interferences to obtain a visible image of an otherwise invisible object. A clear annular gap in the image plane of the condenser is reproduced on the rear image plane of the object lens by the condenser and object lens. A phase shift ring located in the rear image plane of the object lens causes the light passing through this annular gap to shift by a quarter of a wavelength compared to the light diffracted by the object and the part of the rear image plane not covered by the phase shift ring.

The image is generated by the interference effects of the two combined light beams.

Elements of the otherwise invisible samples appear as brightness differences due to the effect of the phase annular gap.

### Installation

To convert the microscope to phase contrast operation, follow the instructions below:

1. Unscrew one of the bright field lenses from the lens turret.

2. Place the desired phase contrast lens into the lens turret.

3. Loosen the fastening screw of the bright field condenser and remove it from the mount. Loosening the locking screw and remove the condenser by pulling it downwards.

4. Insert one of the condensers with the phase contrast diaphragm into the microscope (MML2031 or MML2032)

5. Remove an eyepiece from the tube end and replace it with a centring telescope.

6. Focus the dark pattern of the phase contrast gap of the object lens by adjusting the eyepiece of the centring telescope accordingly.
7. Move the microscope stage with the coarse adjustment knobs to the highest position.
8. Slowly lower the microscope stage until a sharp image of the bright pattern of the condenser annular gap appears.
9. The bright pattern must be adjusted until it is concentric with and overlapped by the darker pattern. Adjust both centring knobs which are located on the front of the ring diaphragm holder of the condenser. The phase annular gap is intentionally made somewhat larger than the image of the condenser diaphragm to ensure a complete overlapping and to prevent light from escaping the phase annular gap. Such a light scattering would affect the contrast.

Non-adjusted phase

Adjusted phase

10. Replace the centring telescope with the previously removed eyepiece.

The microscope can now be used for the examination of samples with the phase contrast method.

*Note!*

The use of an additional Abbe transmitted light condenser whereby a centring diaphragm is mounted for the phase contrast method makes it possible to significantly reduce the time required for converting between bright field and phase contrast method.

## 8. Dark field condenser

The dark field condenser can be mounted to the following microscopes: MML1500, MBL2000 and MBL2000-T. The dark field condenser MML2052 is used for general applications while the dark field condenser MML2053 is especially designed for liquid and blood tests.

### Installation

When changing the microscope for use in dark field operation, take the following steps:

1. Loosen the fastening screw of the bright field condenser and remove it from the mount. Loosening the locking screw and pull the condenser downwards.
2. Now replace the dark field condenser upward into the condenser mount. Adjust the height of the condenser or fasten the locking screw.

The microscope is now ready for use for the dark field operation.

## 9. Photo and video mounting

The trinocular microscope MBL2000-T is equipped with a phototube for mounting photo or video cameras.


The use of the photo adapter MML2042 allows the mounting of a SLR camera to the microscope. The C-mount video adapter MML2047 allows the mounting of a video camera to the microscope.

### Installation

In order to change the microscope to photo or video camera operation, take the following steps:

1. Remove the protective cap from the photo tube.
2. Insert the adapter into the photo tube.
3. When using the photo adapter MML2042, mount the appropriate T2 ring to the adapter. Connect the photo camera to the T2 ring via a bayonet coupling. In case of video adapter MML2047, screw the video camera directly to the adapter.
4. Pull out the light path divider until there is resistance. The light runs through the adapter.
5. The microscope with the photo camera or video camera is now put in operation. In case of a 10x or 40x magnification, a sample is placed below the object lens and the image in the eyepieces is focused using the fine adjustment knob. Now focus the image on the photo camera or video camera. After slightly unscrewing the locking screw of the adapter (S), you can slightly raise or lower the adapter.

## 10. Maintenance

Adjust the position for a sharp image on the camera or the video monitor. This position is fixed by fastening the locking screw on the adapter.

In case of the photo adapter MML2042, it may be necessary to also loosen the stud screw (X) and to change the distance in the lower adapter area to achieve a perfect sharpness.

After these adjustments, the microscope with the photo or video device is ready for use.

## 10. Maintenance

Just like all precision instruments, your microscope should be handled with care. Protect your microscope from dust, condensation, vibration and direct exposure to sunlight.

Keep the optical system very clean. Remove dust with a soft brush or clean air. Remove finger prints with a lint-free cloth soaked in an ether/alcohol mixture (80% ether and 20% alcohol). Special cleaning paper can be purchased in a photo store. Do not use any corrosive solutions for the microscope!

After usage, place the microscope back into the box and protect it with a protective dust cover.

**Important note!**  
Aggressive substances may damage the finish of the plastic casing!

**Important note!**  
In case of a malfunction, the microscope should be repaired by A. KRÜSS Optronic or an authorized repair shop.

## Recovery and Recycling

### Recovery and Recycling

The device can be an important source of raw materials. Please do not dispose of as waste, but collect separately for the recycling and recovery of the contained materials. If disposed improperly, the materials may be damaging to the environment and human health.

The manufacturer of the device, A.KRÜSS Optronic GmbH, collects, uses and recycles the contained raw materials. However this recovery requires your support.

If you decide to dispose of this microscope, please do not try to open it up or to use parts of it in any other way than described in this manual, but return the device to the dealer you purchased it from. The dealer should take the device back free of cost.

The recovery of the raw materials is effected with respect to the European guideline 2002/96/EC and any other applicable guideline.

## Warranty conditions

A. KRÜSS Optronic guarantees the materials and workmanship of KRÜSS microscopes for a period of 24 months from the date of shipping. A. KRÜSS Optronic will repair or replace defective devices within this period, given that they fall under terms of the guarantee. The device must be sent back to A. KRÜSS Optronic for warranty repairs or service. Shipping from the customer is at the expense of A. KRÜSS Optronic for warranty repairs, in any other case at the customer's expense.

A. KRÜSS Optronic guarantees that the hardware specified by A. KRÜSS Optronic for use with this device will function without error if used according to our manufacturer guidelines.

A. KRÜSS Optronic does not guarantee error-free and uninterrupted operation of the device or the accuracy of this instruction manual. A. KRÜSS Optronic is not liable for consequential damage.

### Warranty limitations:

This warranty does not cover errors and damage arising due to improper handling, through modification, misuse, operation above and beyond the specified environmental conditions or through unauthorised maintenance.

Further claims will not be accepted or recognised. A. KRÜSS Optronic expressly provides no guarantee of the workability or economic use in specific application cases.

A. KRÜSS Optronic reserves the right to make changes to this instruction manual and the technical data of the device in question at any point in time.

This digital refractometer may only be transported when properly packed in the complete original packaging including the moulded plastic protectors. Request replacement packaging from your supplier if necessary.

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A. KRÜSS Optronic GmbH Alsterdorfer Straße 276–278 22297 Hamburg | Germany