



Accessories

MA9001 pH 1.68 buffer solution, 230 mL bottle MA9004 pH 4.01 buffer solution, 230 mL bottle MA9006 pH 6.86 buffer solution, 230 mL bottle MA9007 pH 7.01 buffer solution, 230 mL bottle MA9009 pH 9.18 buffer solution, 230 mL bottle pH 10.01 buffer solution, 230 mL bottle MA9010 Refilling Electrolyte Solution 3.5M KCI MA9011 for pH/ORP electrodes, 230 mL Refilling Electrolyte Solution 1M KNO3, MA9012

230 mL, food applications MA9015 Electrode storage solution,

230 mL bottle

MA9016 Electrode cleaning solution, 230 mL bottle MA9112 pH 12.45 buffer solution. . 230 mL bottle 12 VDC Adapter, 220 V MA9310

MA9311 12 VDC Adapter, 110 V MA9315 Electrode Holder

MA917B/1 Glass body, double junction refillable pH electrode

MA831R Temperature probe

MW150 MAX

pH/ORP/Temperature Laboratory Bench Meter

MW150 MAX is an advanced pH/ORP/Temp microprocessor-based bench meter. It is ideal for students and technicians who need fast and reliable measurements.

This meter is provided with a series of new diagnostic features which add an entirely new dimension to the measurement of pH, by allowing the user to dramatically improve the reliability of the measurement:

- Up to 3-point automatic calibration with 7 standard calibration buffers
- Automatic or manual temperature compensation
- Built-in rechargeable battery with 8 hours battery life
- · Auto-off feature to preserve battery energy
- · Battery charger with battery monitor
- · Dedicated GLP key
- Alphanumeric LCD displayed messages for user friendly, intuitive information/warning/ error messages
- Internal clock and date to keep track of different time-dependent functions (calibration, timestamp, calibration time out)
- dedicated GLP key

Glass Electrode & Temperature

Choose from our wide selection of pH and ORP electrodes at pages 8 and



Custom dual level LCD

Large and easyto-read Custom dual level LCD with simultaneous readings and with user-friendly



Ordering Information

MW150 MAX is supplied complete with:

- MA917B/1 Double junction refillable pH electrode
- MA831R Temperature Probe
- MA9315 Electrode Holder
- M10004 pH 4.01 Sachet Buffer Solution • M10007 pH 7.01 Sachet Buffer Solution
- M10010 pH 10.01 Sachet Buffer Solution
- M10016 Sachet Electrode Cleaning Solution

rehno.com

- MA9310 12 VDC Adapter
- · Graduated pipette
- Instruction manual

